


PRIMEDIfem®

PRODUCT CATALOG



primed®



Gentle and effektive

PRIMEDIfem® Intrauterine Contraceptive Devices

PRIMEDIfem® IUDs are high-quality products suitable for all women who desire a simple and hormone free method of contraception which remains effective over several years.

HORMON-FREE:

The hormone-free PRIMEDIfem® IUDs avoid the typical risks of long-term hormonal-contraception, such as cellular changes in breast and ovaries, vascular diseases (thromboembolism) and cardiovascular diseases.¹

SAFE AND SECURE:

Contraception by means of an IUD is one of the most reliable methods around. Worldwide more than 100 million women have a small plastic device with a copper coil (pessary with copper coil) inserted in their uterus to avoid unwanted pregnancy.

Not only does the IUD itself act as a contraceptive, the copper also boosts this effect. The copper coil evenly releases minute quantities of copper into the uterine mucosa with the result that the mobility of the sperm is restricted to the point that they cannot enter the Fallopian tubes to fertilise the ovum.

The contraceptive effect begins immediately after insertion. Fertility is restored immediately after removing the IUD.


PRIMEDIfem® IUDs feature a copper wire of 99,99 % purity. Undesireble allergic reactions can therefore be excluded.

The return threads made from monofilament, non-absorbable synthetics are tear-proof and yet soft.

1. T. Soini, R. Hurskainen, S. Grénman, J. Mäenpää, J. Paavonen, H. Joensuu, E. Pukkala: Levonorgestrel-releasing intrauterine system and the risk of breast cancer: A nationwide cohort study. In: Acta oncologica. Band 55, Nummer 2, 2016, S. 188–192; Ärzteblatt über psychiatrische Nebenwirkungen der Hormonspirale, 2009.


DIFFERENT DESIGNS:

PRIMEDIfem® IUDs are available in many different shapes and versions.



U-shaped


- Well-established and trusted model
- Secure fit thanks to the curved flexible side arms
- Available with additional goldclips



T-shaped

- Flexible, lily-shaped polyethylene frame
- Quick-load insertion system with applicator
- Easy removal
- Available in two different sizes, normal and mini (for small uterine cavities)
- Available with additional goldclips

PRIMEDIfem® CU 375 Intrauterine Contraceptive Device



PRIMEDIfem® CU 375

- Copper surface: 375 mm²
- In-situ-time up to 5 years
- Secure fit in the uterus
- Sterile

PU/SU


Article number

1/10

20146

PU = packaging unit, SU = Shipping unit

PRIMEDIfem® CU 375 gold Intrauterine Contraceptive Device



PRIMEDIfem® CU 375 gold

- Copper surface: 375 mm²
- In-situ-time up to 5 years
- Secure fit in the uterus
- With additional goldclips
- Sterile


PU/SU

Article number

1/10

20138

PU = packaging unit, SU = Shipping unit

PRIMEDIfem[®]-T CU 380 Intrauterine Contraceptive DevicePRIMEDIfem[®]-T CU 380

- Copper surface: 380 mm²
- In-situ-time up to 5 years
- Secure fit in the uterus
- Sterile

PU/SU

Article number

1/10

20145

PU = packaging unit, SU = Shipping unit

PRIMEDIfem[®]-T mini CU 380 Intrauterine Contraceptive DevicePRIMEDIfem[®]-T mini CU 380

- Copper surface: 380 mm²
- In-situ-time up to 5 years
- Secure fit in the uterus
- Suitable for women with a small uterine cavity
- Sterile

PU/SU

Article number

1/10

20147

PU = packaging unit, SU = Shipping unit

PRIMEDIfem[®]-T CU 380 gold Intrauterine Contraceptive DevicePRIMEDIfem[®]-T CU 380 gold

- Copper surface: 380 mm²
- In-situ-time up to 5 years
- Secure fit in the uterus
- With additional goldclips
- Sterile

PU/SU


Article number

1/10

20214

PU = packaging unit, SU = Shipping unit

PRIMEDIfem[®]-T mini CU 380 gold Intrauterine Contraceptive Device



PRIMEDIfem[®]-T mini CU 380 gold

- Copper surface: 380 mm²
- In-situ-time up to 5 years
- Secure fit in the uterus
- Suitable for women with a small uterine cavity
- With additional goldclips
- Sterile

PU/SU	Article number
1/10	20216

PU = packaging unit, SU = Shipping unit

Manufacturer

Primed[®] Halberstadt
Medizintechnik GmbH
Strasse des 20. Juli 1
D-38820 Halberstadt
Phone: 0049-(0)3941-668-6
Fax: 0049-(0)3941-24565
primed@primed-halberstadt.de
www.primed-halberstadt.de

Distributor